


URETERO-RENSCOPE



MANUFACTURED BY

APPLE LIFE SCIENCES

W-19/15 Western Avenue, Sainik Farms, New Delhi - 110 062

Ph. : 011-29551678 +91 9810264605 +91 9711484602

info@applelifesciences.com

www.applelifesciences.com

URETERO-RENSCOPIES AND ACCESSORIES , APPLE MADE IN GERMANY

CAT. NO.	DESCRIPTION
25-01-4425	HIGH DEFINATION, COMPACT FIBRE URETERO-RENSCOPE, AUTOCLAVABLE, DIRECTION OF VIEW 0°, OPTICAL FIBRE BUNDLE 17000 PIXELS. DISTAL TIP 4.5 FR., OUTER TUBE Ø 6.5 FR. – 8.1 FR., WORKING LENGTH 425 MM. LARGE IRRIGATION AND INSTRUMENT CHANNEL 4 FR. FOR USE WITH 1X3 FR. AUXILIARY INSTRUMENTS IT HAS THREE PLASTIC MAINTENANCE FREE STOP COCKS , LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF IRRIGATION FLUID
25-215-3600	STONE GRASPING FORCEPS , SEMI-RIGID, FOR URETERO-RENSCOPE, DIA, 3 FR., WORKING LENGTH 600 MM.
25-220-3600	BIOPSY FORCEPS, SEMI-RIGID, OVAL CUP SHAPED JAWS, FOR URETERO-RENSCOPE, DIA, 3 FR., WORKING LENGTH 600 MM.
10-180-3600	COAGULATION BALL ELECTRODE FOR URETEROSCOPE DIA 3 FR. WORKING LENGTH 600 MM.
10-170-3600	COAGULATION NEEDLE ELECTRODE FOR URETEROSCOPE DIA 3 FR. WORKING LENGTH 600 MM
25-01-4310	HIGH DEFINATION, COMPACT FIBRE URETERO-RENSCOPE SHORT, AUTOCLAVABLE, DIRECTION OF VIEW 0°. OPTICAL FIBRE BUNDLE 17000 PIXELS. DISTAL TIP 4.5 FR., OUTER TUBE Ø 6.5 FR. – 8.1 FR., WORKING LENGTH 310 MM. LARGE IRRIGATION AND INSTRUMENT CHANNEL 4 FR. FOR USE WITH 1X3 FR. AUXILIARY INSTRUMENTS IT HAS THREE PLASTIC MAINTENANCE FREE STOP COCKS , LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF IRRIGATION FLUID
25-215-3425	STONE GRASPING FORCEPS, SEMI-RIGID, FOR SHORT URETERO-RENSCOPE, DIA, 3 FR., WORKING LENGTH 425 MM
25-01-6425	HIGH DEFINATION, COMPACT FIBRE URETERO-RENSCOPE, AUTOCLAVABLE, DIRECTION OF VIEW 0°, OPTICAL FIBRE BUNDLE 17000 PIXELS. DISTAL TIP 6 FR., OUTER TUBE Ø 7.5 FR. - 9 FR., WORKING LENGTH 425 MM. LARGE IRRIGATION AND INSTRUMENT CHANNEL 5.1 FR. FOR USE WITH 1X4 FR. OR 2 X 2.4 FR. AUXILIARY INSTRUMENTS IT HAS THREE PLASTIC MAINTENANCE FREE STOP COCKS , LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF IRRIGATION FLUID
25-204-4600	STONE GRASPING FORCEPS, SEMI-RIGID, ALLIGATOR JAWS, FOR URETERO-RENSCOPE, DIA, 4 FR., WORKING LENGTH 600 MM
25-216-4600R	THREE PRONGED GRASPING FORCEPS, WITH RING HANDLE, FOR URETERO-RENSCOPE, DIA, 4 FR., WORKING LENGTH 600 MM.
25-218-4600	BIOPSY FORCEPS, SEMI-RIGID, OVAL CUP SHAPED JAWS. FOR URETERO-RENSCOPE, DIA, 4 FR., WORKING LENGTH 600 MM.
25-224-4600	SCISSOR, SINGLE ACTION JAWS, SEMI-RIGID, FOR URETERO-RENSCOPE DIA, 4 FR., WORKING LENGTH 600 MM.
10-180-4600	COAGULATION BALL ELECTRODE FOR URETEROSCOPE DIA 4 FR. WORKING LENGTH 600 MM.
10-170-4600	COAGULATION NEEDLE ELECTRODE FOR URETEROSCOPE DIA 4 FR. WORKING LENGTH 600 MM
25-01-6310	HIGH DEFINATION, COMPACT FIBRE URETERO-RENSCOPE SHORT, AUTOCLAVABLE, DIRECTION OF VIEW 0°. OPTICAL FIBRE BUNDLE 17000 PIXELS. DISTAL TIP 6 FR., OUTER TUBE Ø 7.5 FR. - 9 FR., WORKING LENGTH 310 MM. LARGE IRRIGATION AND INSTRUMENT CHANNEL 5.1 FR. FOR USE WITH 1X4 FR. OR 2 X 2.4 FR. AUXILIARY INSTRUMENTS IT HAS THREE PLASTIC MAINTENANCE FREE STOP COCKS , LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF IRRIGATION FLUID
25-204-4425	STONE GRASPING FORCEPS, SEMI-RIGID, ALLIGATOR JAWS, FOR SHORT URETERO-RENSCOPE, DIA, 4 FR., WORKING LENGTH 425 MM
25-01-8425	HIGH DEFINATION, COMPACT FIBRE URETERO-RENSCOPE, AUTOCLAVABLE, DIRECTION OF VIEW 0°, OPTICAL FIBRE BUNDLE 40000 PIXELS. DISTAL TIP 7.5 FR., OUTER TUBE Ø 9 FR. - 12 FR., WORKING LENGTH 425 MM. LARGE IRRIGATION AND INSTRUMENT CHANNEL 6.8 FR. FOR USE WITH 1X5 FR. OR 2 X 2.4 FR. AUXILIARY INSTRUMENTS IT HAS THREE PLASTIC MAINTENANCE FREE STOP COCKS , LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF IRRIGATION FLUID
	25-204-5600 STONE GRASPING FORCEPS, SEMI-RIGID, KING SIZE ALLIGATOR JAWS, FOR URETERO-RENSCOPE, DIA, 5 FR., WORKING LENGTH 600 MM
25-212-5600	STONE GRASPING FORCEPS, SEMI-RIGID, ALLIGATOR JAWS. FOR URETERO-RENSCOPE, DIA, 5 FR., WORKING LENGTH 600 MM.
25-214-5600	THREE PRONGED GRASPING FORCEPS, FOR URETERO-RENSCOPE, DIA, 5 FR., WORKING LENGTH 600 MM.
25-216-5600R	THREE PRONGED GRASPING FORCEPS, WITH RING HANDLE, FOR URETERO-RENSCOPE, DIA, 5 FR., WORKING LENGTH 600 MM.
25-220-5600	BIOPSY FORCEPS, SEMI-RIGID, OVAL CUP SHAPED JAWS. FOR URETERO-RENSCOPE DIA, 5 FR., WORKING LENGTH 600 MM.
25-224-5600	SCISSOR, SEMI-RIGID, SINGLE ACTION JAWS. FOR URETERO-RENSCOPE, DIA, 5 FR., WORKING LENGTH 600 MM.
10-180-5600	COAGULATION BALL ELECTRODE FOR URETEROSCOPE DIA 5 FR. WORKING LENGTH 600 MM.
10-170-5600	COAGULATION NEEDLE ELECTRODE FOR URETEROSCOPE DIA 5 FR. WORKING LENGTH 600 MM
25-01-8310	HIGH DEFINATION, COMPACT FIBRE URETERO-RENSCOPE SHORT, AUTOCLAVABLE, DIRECTION OF VIEW 0°. OPTICAL FIBRE BUNDLE 40000 PIXELS. DISTAL TIP 7.5 FR., OUTER TUBE Ø 9 FR. - 12 FR., WORKING LENGTH 310 MM. LARGE IRRIGATION AND INSTRUMENT CHANNEL 6.8 FR. FOR USE WITH 1X5 FR. OR 2 X 3 FR. AUXILIARY INSTRUMENTS IT HAS THREE PLASTIC MAINTENANCE FREE STOP COCKS , LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF IRRIGATION FLUID
25-204-5425	STONE GRASPING FORCEPS, SEMI-RIGID, KING SIZE ALLIGATOR JAWS, FOR SHORT URETERO-RENSCOPE, DIA, 5 FR., WORKING LENGTH 425 MM
25-202-5425	STONE GRASPING FORCEPS, SEMI-RIGID, FOR SHORT URETERO-RENSCOPE, DIA, 5 FR., WORKING LENGTH 425 MM