



MINI NEPHROSCOPE 12 FR



MANUFACTURED BY

APPLE LIFE SCIENCES

W-19/15 Western Avenue, Sainik Farms, New Delhi - 110 062

Ph. : 011-29551678 +91 9810264605 +91 9711484602

info@applelifesciences.com

www.applelifesciences.com

MINI NEPHROSCOPE 12 FR., APPLE, MADE IN GERMANY.

CAT. NO.	DESCRIPTION
25-12-1222	COMPACT MINI NEPHROSCOPE, AUTOCLAVABLE, DIRECTION OF VIEW 0°, Ø 12 FR., OPTICAL FIBRE BUNDLE 40000 PIXELS. WITH OBLIQUE EYE PIECE. WORKING LENGTH 220 MM, STRIGHT INSTRUMENT CHANNEL 6.6 FR. FOR USE WITH 1 X5 FR INSTRUMENT. IT HAS TWO PLASTIC MAINTENANCE FREE STOP COCKS, LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF THE IRRIGATION FLUID.
25-15-015	OPERATING SHEATH FOR MINI NEPHROSCOPE FOR CONTINUOUS IRRIGATION AND ASPIRATION, DIAMETER 15 FR., WORKING LENGTH 150 MM.
25-22-015	ONE STEP DILATOR WITH CENTRAL CHANNEL FOR GUIDE WIRE TO BE USED WITH 15 FR. OPERATION SHEATH FOR MINI NEPHROSCOPE, WORKING LENGTH 300 MM.
25-15-016	OPERATING SHEATH FOR MINI NEPHROSCOPE FOR CONTINUOUS IRRIGATION AND ASPIRATION, DIAMETER 16.5 FR., WORKING LENGTH 150 MM.
25-22-016	ONE STEP DILATOR WITH CENTRAL CHANNEL FOR GUIDE WIRE TO BE USED WITH 16.5 FR. OPERATION SHEATH FOR MINI NEPHROSCOPE, WORKING LENGTH 300 MM.
25-15-021	OPERATING SHEATH FOR MINI NEPHROSCOPE FOR CONTINUOUS IRRIGATION AND ASPIRATION, DIAMETER 21 FR., WORKING LENGTH 150 MM.
25-22-021	ONE STEP DILATOR WITH CENTRAL CHANNEL FOR GUIDE WIRE TO BE USED WITH 21 FR. OPERATION SHEATH FOR MINI NEPHROSCOPE, WORKING LENGTH 300 MM.
25-14-5410	STONE GRASPING FORCEPS, SEMI-RIGID, ALLIGATOR JAWS FOR MINI NEPHROSCOPE Ø 5 FR., WORKING LENGTH 410 MM.
25-16-5410	BIOPSY FORCEPS, OVAL CUP SHAPED JAWS, SEMI-RIGID, Ø 5 FR., WORKING LENGTH 410 MM
25-18-5410	MICRO SCISSOR, SINGLE ACTION JAWS, SEMI-RIGID, Ø 5 FR., WORKING LENGTH 410 MM
25-26-011	THREE PRONGED, STONE GRASPING FORCEPS, U-SHAPED SPRING HANDLE, Ø 5 FR., WORKING LENGTH 360 MM