



NEPHROSCOPE 18.8 FR



MANUFACTURED BY

APPLE LIFE SCIENCES

W-19/15 Western Avenue, Sainik Farms, New Delhi - 110 062

Ph. : 011-29551678 +91 9810264605 +91 9711484602

info@applelifesciences.com

www.applelifesciences.com

COMPACT NEPHROSCOPE 18.8 FR. APPLE, MADE IN GERMANY	
CAT. NO.	DESCRIPTION
25-10-1822	COMPACT NEPHROSCOPE, PARALLEL EYE PIECE, AUTOCLAVABLE, Ø 18.8 FR., DIRECTION OF VIEW 15°, INSTRUMENT CHANNEL 4.3 MM. FIELD OF VIEW 80°± 5%, WORKING LENGTH 220 MM.
25-13-1822	COMPACT NEPHROSCOPE, OBLIQUE EYE PIECE, AUTOCLAVABLE, Ø 18.8 FR., DIRECTION OF VIEW 15°, INSTRUMENT CHANNEL 4.3 MM. FIELD OF VIEW 80°± 5%, WORKING LENGTH 220 MM. INCLUDES IRRIGATION ADAPTOR WITH STOPCOCK (25-15-000)
25-15-000	IRRIGATION ADAPTOR FOR COMPACT NEPHROSCOPE WITH STOP COCK.
25-15-001	SILICON VALVE FOR IRRIGATION ADAPTOR
25-15-024	OPERATING SHEATH WITH OBTURATOR AND SWIVEL IRRIGATION CONNECTOR, CONTINUOUS IRRIGATION SIZE 24 FR.
25-24-3538	STONE GRASPING FORCEPS WITH ALLIGATOR JAWS, Ø 3.5 MM. WORKING LENGTH 380 MM
25-26-3538	PEANUT STONE GRASPING FORCEPS, DOUBLE ACTION JAWS WITH FENESTRATION, Ø 3.5 MM. WORKING LENGTH 380 MM
25-28-3538	STONE GRASPING FORCEPS WITH FENESTRATED JAWS, U-SHAPED SPRING HANDLE, Ø 3.5 MM, WORKING LENGTH 380 MM.
25-30-3538U	THREE PRONGED, STONE GRASPING FORCEPS, U-SHAPED SPRING HANDLE, Ø 3.5 MM, WORKING LENGTH 380 MM.