

PAEDIATRIC COMPACT CYSTO-URETHEROSCOPE



MANUFACTURED BY

APPLE LIFE SCIENCES

W-19/15 Western Avenue, Sainik Farms, New Delhi - 110 062

Ph. : 011-29551678 +91 9810264605 +91 9711484602

info@applelifesciences.com

www.applelifesciences.com

PAEDIATRIC CYSTO-URETHROSCOPE APPLE MADE IN GERMANY

CAT. NO.	DESCRIPTION
22-01-4150	COMPACT FIBRE CYSTO-URETHROSCOPE HIGH DEFINITION, AUTOCLAVABLE, DIRECTION OF VIEW 0 DEGREE, OPTICAL FIBRE BUNDLE 17000 PIXELS. DISTAL TIP 4.5/6.5 FR., WORKING LENGTH 150 MM, LARGE IRRIGATION AND INSTRUMENT CHANNEL 4 FR FOR USE WITH 1X3 FR. INSTRUMENTS, IT HAS THREE PLASTIC MAINTENANCE FREE STOP COCKS, LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF THE IRRIGATION FLUID.
10-150-003	GRASPING FORCEPS, FLEXIBLE WITH ALLIGATOR JAWS, Ø 3 FR., WORKING LENGTH 350 MM
10-160-003	BIOPSY FORCEPS, FLEXIBLE, OVAL CUP JAWS, Ø 3 FR., WORKING LENGTH 350 MM
10-168-003	SCISSOR, FLEXIBLE, SINGLE ACTION JAWS, Ø 3 FR., WORKING LENGTH 350 MM
10-170-3350	BUGBY NEEDLE ELECTRODE, MONO-POLAR, 3 FR., 350 MM LONG
10-180-3350	BUGBY BALL ELECTRODE, MONO-POLAR, 3 FR., 350 MM LONG
10-410-03	HF CABLE MONO-POLAR FOR ABOVE
22-01-6150	COMPACT FIBRE CYSTO-URETHROSCOPE HIGH DEFINITION, AUTOCLAVABLE, DIRECTION OF VIEW 0 DEGREE, OPTICAL FIBRE BUNDLE 17000 PIXELS. DISTAL TIP 6/7.5-9 FR., WORKING LENGTH 150 MM LARGE IRRIGATION AND INSTRUMENT CHANNEL 5.1 FR., TO BE USED FOR 1X4 OR 2X2.4 FR. INSTRUMENTS. IT HAS THREE PLASTIC MAINTENANCE FREE STOP COCKS, LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF THE IRRIGATION FLUID.
25-14-4400	STONE GRASPING FORCEPS, SEMI-RIGID, ALLIGATOR JAWS FOR COMPACT CYSTO-URETHROSCOPE 6/7.5 FR. Ø 4 FR., WORKING LENGTH 400 MM.
10-162-004	BIOPSY FORCEPS, SEMI-RIGID, OVAL CUP JAWS, FOR COMPACT CYSTO-URETHROSCOPE 6-7.5 – 9 FR. Ø 4 FR. WORKING LENGTH 400 MM
10-170-4350	BUGBY NEEDLE ELECTRODE, MONO-POLAR, 4 FR., 350 MM LONG
10-180-4350	BUGBY BALL ELECTRODE, MONO-POLAR, 4 FR., 350 MM LONG
10-410-03	HF CABLE MONO-POLAR FOR ABOVE
22-01-8150	COMPACT FIBRE CYSTO-URETHROSCOPE HIGH DEFINITION, AUTOCLAVABLE, DIRECTION OF VIEW 0 DEGREE, OPTICAL FIBRE BUNDLE 40000 PIXELS. DISTAL TIP 7.5/9.8 FR., WORKING LENGTH 150 MM, LARGE IRRIGATION AND INSTRUMENT CHANNEL 6.8 FR., TO BE USED FOR 1X5 FR. OR 2 X3 FR. INSTRUMENT. IT HAS THREE PLASTIC MAINTENANCE FREE STOP COCKS, LUER CONNECTOR AND STRAIGHT WORKING CHANNEL WITH SILICON VALVE TO STOP BACK FLOW OF THE IRRIGATION FLUID.
25-14-5400	STONE GRASPING FORCEPS, SEMI-RIGID, ALLIGATOR JAWS FOR COMPACT CYSTO-URETHROSCOPE 7.5-9 - 12 FR. Ø 5 FR. WORKING LENGTH 400 MM.
10-162-005	BIOPSY FORCEPS, SEMI-RIGID, OVAL CUP JAWS, FOR COMPACT CYSTO-URETHROSCOPE 7.5-9 - 12 FR. Ø 5 FR. WORKING LENGTH 400 MM
10-170-5350	BUGBY NEEDLE ELECTRODE, MONO-POLAR, 5 FR., 350 MM LONG
10-180-5350	BUGBY BALL ELECTRODE, MONO-POLAR, 5 FR., 350 MM LONG
10-410-03	HF CABLE MONO-POLAR FOR ABOVE