



# PAEDIATRIC RESECTOSCOPE



**MANUFACTURED BY**

**APPLE LIFE SCIENCES**

W-19/15 Western Avenue, Sainik Farms, New Delhi - 110 062

Ph. : 011-29551678 +91 9810264605 +91 9711484602

info@applelifesciences.com

www.applelifesciences.com

<b>PAEDIATRIC RESECTOSCOPE, APPLE MADE IN GERMANY</b>	
<b>CAT. NO.</b>	<b>DESCRIPTION</b>
22-10-190	FIBER OPTIC STRAIGHT FORWARD TELESCOPE 0°, Ø 1.9 MM, 187 MM LONG WITH LATEST IMAGE GUIDED TECHNOLOGY USING SPECIAL SILICON FIBERS FOR OPTIMUM IMAGE QUALITY AND LONGER LIFE
22-10-193	FIBER OPTIC STRAIGHT FORWARD TELESCOPE 30°, Ø 1.9 MM, 187 MM LONG WITH LATEST IMAGE GUIDED TECHNOLOGY USING SPECIAL SILICON FIBERS FOR OPTIMUM IMAGE QUALITY AND LONGER LIFE
22-10-009	RESECTOSCOPE SHEATH 9 FR. WITH FIXED IRRIGATION TAP, DISTAL TIP INSULATED AND OBTURATOR (22-10-019)
22-12-009	RESECTOSCOPE WORKING ELEMENT, PASSIVE
22-15-101	PAEDIATRIC CUTTING LOUPE FOR 9 FR. RESECTOSCOPE
22-15-102	PAEDIATRIC HOOK ELECTRODE FOR 9 FR. RESECTOSCOPE
22-15-103	PAEDIATRIC COAGULATION ELECTRODE FOR 9 FR. RESECTOSCOPE
10-410-03	HF CABLE MONO-POLAR FOR ABOVE