

# RESECTOSCOPE CONTINUOUS IRRIGATION



**MANUFACTURED BY**

**APPLE LIFE SCIENCES**

W-19/15 Western Avenue, Sainik Farms, New Delhi - 110 062

Ph. : 011-29551678 +91 9810264605 +91 9711484602

info@applelifesciences.com

www.applelifesciences.com

**MONO-POLAR RESECTOSCOPE CONTINUOUS IRRIGATION 26/24 FR.  
APPLE, MADE IN GERMANY**

<b>CAT. NO.</b>	<b>DESCRIPTION</b>
10-100-030	HD TELESCOPE FORWARD OBLIQUE 30°, WIDE ANGLE, Ø 4 MM, 303 MM LONG, AUTOCLAVABLE, COLOUR CODE: RED
10-300-026	RESECTOSCOPE OUTER SHEATH, CONTINUOUS IRRIGATION, DOUBLE SHEATH SYSTEM, SIZE 26 FR.
10-300-024	RESECTOSCOPE INNER SHEATH WITH OBLIQUE CERAMIC TIP, ROTATABLE 360°, SIZE 24 FR.
10-300-124	STANDARD OBTURATOR FOR RESECTOSCOPE SHEATH 24 FR.
10-332-024	DEFLECTING OBTURATOR FOR RESECTOSCOPE SHEATH 24 FR.
10-334-024	VISUAL OBTURATOR FOR RESECTOSCOPE SHEATH 24 FR.
10-350-024	WORKING ELEMENT FOR RESECTOSCOPE, ACTIVE (CUTTING BY MEANS OF FINGER GRIP AND THUMB IS STATIONARY)
10-352-024	WORKING ELEMENT FOR RESECTOSCOPE, PASSIVE (CUTTING BY MEANS OF SPRING ACTION, THUMB IS MOVING AND FINGER IS STATIONARY)
10-400-100	CUTTING LOUPE ANGLED Ø 0.35 MM FOR RESECTOSCOPE 26/24 FR.
10-400-101	CUTTING LOUPE STRAIGHT Ø 0.35 MM FOR RESECTOSCOPE 26/24 FR.
10-400-102	COAGULATION ELECTRODE, CONICAL FOR RESECTOSCOPE
10-400-103	COAGULATION BALL ELECTRODE Ø 3 MM FOR RESECTOSCOPE
10-400-104	COAGULATION ROLLER BALL ELECTRODE Ø 5 MM FOR RESECTOSCOPE
10-400-105	COLLIN'S KNIFE FOR RESECTOSCOPE
10-430-000	ELLIKEVACUATOR WITH STANDARD CONNECTION
10-410-03	MONO-POLAR HF CONNECTING CABLE FOR RESECTOSCOPE
	<b>BI-POLAR RESECTOSCOPE CONTINUOUS IRRIGATION 26/24 FR. APPLE, MADE IN GERMANY</b>
10-100-030	HD TELESCOPE FORWARD OBLIQUE 30°, WIDE ANGLE, Ø 4 MM, 303 MM LONG, AUTOCLAVABLE, COLOUR CODE: RED
10-300-026	RESECTOSCOPE OUTER SHEATH, CONTINUOUS IRRIGATION, SIZE 26 FR.
10-300-024	RESECTOSCOPE INNER SHEATH WITH OBLIQUE CERAMIC TIP, ROTATABLE 360°, SIZE 24 FR.
10-330-024	STANDARD OBTURATOR FOR RESECTOSCOPE SHEATH 24 FR.
10-332-024	DEFLECTING OBTURATOR FOR RESECTOSCOPE SHEATH 24 FR.
10-334-024	VISUAL OBTURATOR FOR RESECTOSCOPE SHEATH 24 FR.
10-354-024	BI-POLAR WORKING ELEMENT, ACTIVE
10-356-024	BI-POLAR WORKING ELEMENT, PASSIVE
10-410-100	BI-POLAR CUTTING LOUPES, ANGLED, Ø 0.35 MM FOR RESECTOSCOPE 26/24 FR.
10-410-103	BI-POLAR COAGULATION ELECTRODE FOR RESECTOSCOPE
10-410-105	BI-POLAR COLLIN'S KNIFE FOR RESECTOSCOPE
10-412-03	BI-POLAR CABLES FOR APPLE RESECTOSCOPE